




SUCTION FITTINGS SYSTEM






Index

Place selector valves.....	3-4
Female connectors for hose	5
Male connectors for hose.....	6
Female-Female straight connectors	7
Female-Male straight connectors	8-9
Special connectors.....	10
T-junctions.....	11
Secretion filter.....	12-26
Disposable filter.....	27
O-ring.....	28
Suction handpieces.....	29-30
Aluminium hose adapter.....	31-32
Universal pvc hose adapter.....	33-34
Suction handpieces and parts.....	35
Handpiece holders.....	36
Complete suction hoses.....	37-38
Suction hoses.....	39-41
Suction release valve.....	42

Place selector valves



Place selector valves





DE5.19D

24V AC/DC
Connectors 25 mm





DE5.19DR

24V AC/DC
Connectors 25 mm







F connectors for hose




TM72006E




TM72009E




TM72010E




TM72012E




TM72013E




TM72014E




TM72016E




TM72018E




TM72019E




TM72020E



TM72021E



TM72022E




TM72025E



TM72027E



M connectors for hose




TM71006




TM71008




TM71010




TM71012




TM71013



TM71014



TM71016



TM71018



TM71019



TM71020



TM71021



TM71022





TM71025




TM71027





F-F straight connectors





TM70011E




TM70014E



TM70028E





TM70028EB







F-M straight connectors




TM70017E





TM70015E




TM70043E




TM70018E




TM70018EB



F-M straight connectors



4001TM22



Ø22



Ø10.6

4001TM000



Ø11.8

4001TM001




Ø18


4001TM002




Special connectors




TM70004E




TM70005E




TM101010



TM70001
not included





TM80010E




TM70051E




T-Junctions




TM80006E




TM70008E



TM70003E




DE5.31.1



TM70006E


Secretion filter

Page 5




DE1002X.15

CLOSE




152

70



Ø10.6

4001TM000




Ø11.8

4001TM001



Ø18

4001TM002




72

Ø37


150

DE1002X.15.1 replacement filter

Secretion filter



Page 5




TM702.055

CLOSE



Ø10.6

4001TM000




Ø11.8

4001TM001





Ø18

4001TM002



152


70








TM725.044 replacement filter





Connection example







Secretion filter



Secretion filter




flexible
silicone
plug









Connection example







Secretion filter

Page 5





Secretion filter




Page 5




Secretion filter



CLOSE







Ø6.5 TCM065D
 Ø8.5 TCM085D
 Ø9.5 TCM095D
 Ø11 TCM110D
 Ø12 TCM120D
 Ø15 TCM150D
 Ø18 TCM180D




Secretion filter



Page 5




Secretion filter




TMY70270E (YELLOW)
TMG70270E (GRAY)





Secretion filter



TMY70272E (YELLOW)
TMG70272E (GRAY)





Ø6.5 TCM065D Ø8.5 TCM085D Ø9.5 TCM095D Ø11 TCM110D Ø12 TCM120D Ø15 TCM150D Ø18 TCM180D








Connection example




Secretion filter

DE5.15.1












Secretion filter





Disposable filters



O-ring



OR1602/G



OR1602/P



OR3230




flexible
silicone
plug

No o-ring




OR2118



OR3352

OR1815



OR2118

OR2802


Suction handpieces





Suction handpieces

040T515.1



040T515
Ø5.5/7



flexible silicone

135°C
SSS

040T410P
Ø6.5



flexible silicone

135°C
SSS

040T415P
Ø6.5



135°C
SSS

040T4102P
2X Ø6.5



135°C
SSS

040T520P

Ø11



flexible silicone

135°C
SSS


040T515
Ø5.5/7



flexible silicone

135°C
SSS

Hose connector included in the VC1602/P VC1602PA



Ø11.1

2T5238



VC1602/P



VC1602PA


Ø5.5/7



VC1604P

flexible silicone

135°C
SSS



2T5239



VC1602/G



VC1602GA



VC1604G


135°C
SSS

135°C
SSS

135°C
SSS

135°C
SSS

Hose connector included in the VC1602/G VC1602GA



Ø17.5



Aluminium hose adapter



Ø16.5





Aluminium hose adapter



Ø22.5



TCM0013
Isodent®/Castellini®



TCM0012
Planmeca®



N.C.
Cattani®



1234251
Sirona®



1234252
Eurodent®



1234255
KaVo®



Ø 17.5 mm

AS17TB10R
10 mt



Ø 17.5 mm

US17TB10R
10 mt



Ø 15 mm

AS15TB10R 10 mt



TB15E.1
Eurodent®




Ø 20 mm

AS20TB10R
10 mt



Universal pvc hose adapter



Universal pvc hose adapter




VC1602/G VC1602GA VC1604G

Ø22.5

TCM110GC



Ø11



Ø11 mm


AS11TB10R

10 mt

TCM120GC



Ø12



Ø13 mm

US11TB10R

10 mt

TCM155GC



Ø15.5



Ø15 mm

AS15TB10R

10 mt

TCM152GD



Ø15.2



Ø17.5 mm


AS17TB10R

10 mt

TCM175GD



Ø17.5




Ø17.5 mm

US17TB10R

10 mt

Suction handpieces and parts



Suction terminal code	Compatibility	Material	Weight (gr)	1 cannula Ø insertion	2 Can be equipped with saliva ejector adaptor ø 6,5 mm	3 Saliva ejector adaptor code	3 Suction terminal lever code	4 O-ring code	5 Ø Tube terminal adaptor	5 Tube terminal adaptor code
VC1602/P VC1602PA*	Cattani®	PVC *PVC	17	11 mm	●	040T515 (Opt)	2T5238	OR1602/P	I 10,8 mm	
VC1604P	Cattani®	Aluminium	29	6,5 mm	●	040T515.1	121T9300	OR1604P (Kit of 4 pc.)	Quick connector Ø16 mm	
VC1604PC	Isodent®/ Castellini®		43						I 10,2 mm	TCM0014
VC1604PE	Eurodent®		38						I 8,4 mm	1234254
VC1604PK	Kavo®		44						I 8,4 mm	TCM0010
VC1604PP	Plameca®		56						I 10,5 mm	TCM0011
VC1604PS	Sirona®		41						I 8,4 mm	1234253
VC1604PD	Dci®		32						I 3/16"	
VC1602/G VC1602GA*	Cattani®		PVC *PVC						27	11-16 mm
VC1604G	Cattani®	Aluminium	61	11 mm	●	040T515 (Opt)	111T7400	OR1604G (Kit of 5 pc.)	Quick connector Ø22 mm	
VC1604GC	Isodent®/ Castellini®		84						I 17,2 mm	TCM0013
VC1604GE	Eurodent®		72						I 14,7 mm	1234252
VC1604GK	Kavo®		89						I 20 mm	1234255
VC1604GP	Plameca®		78						I 17,2 mm	TCM0012
VC1604GS	Sirona®		76						I 14,7 mm	1234251
VC1604GD	Dci®		63						I 1/2"	



Handpiece holders



033T316




033T331



033T333

Complete suction hoses



For Eurodent 8 mm




For Kavo 8 mm



For Sirona 8 mm




For Cattani 11 mm




For Castellini 11 mm



For Planmeca 11 mm



Complete suction hoses




For Eurodent 15 mm



For Kavo 20 mm



For Sirona 15 mm




For Cattani 17 mm



For Castellini 17 mm



For Planmeca 17 mm



Suction hoses ultraflex-d



L = 10 mt


- 15°C to + 70°C

Code	inner d.	outside d.	Wall thickness	vacuum
for 10 mt.	mm	mm	mm	mbar
AS08TB10R	8	12	2.0	250
AS11TB10R	11	15	2.0	240
AS15TB10R	15	19	2.0	220
AS17TB10R	17	21	2.0	200
AS20TB10R	20	25	2.5	180
DE1002X.11	40	46	3.0	150

Suction hoses ultrasoft-d



L = 10 mt



- 5°C to + 60°C

Code	inner d.	outside d.	Wall thickness	vacuum
for 10 mt.	mm	mm	mm	mbar
US11TB10R	13	16	1.5	200
US17TB10R	17	20	1.5	200

Technical pipes



Code	inner d.	vacuum
for 10 mt.	mm	mbar
040T920	15	950
040T950	20	950
040T990	25	950

Suction release valve

