








DENTAL HOSES

| | |
|--|-------------|
|  CEFLA | <u>P. 2</u> |
|  OPEN END TURBINE | <u>P.14</u> |
|  OPEN END MICROMOTOR | <u>P.19</u> |
|  OPEN END SCALER | <u>P.29</u> |
|  KAVO | <u>P.33</u> |
|  PLANMECA | <u>P.44</u> |
|  SIRONA | <u>P.53</u> |

CEFLA DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty










Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N










MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--------------------------|---|---|--|---|------------|
| MICROMOTOR MC2 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC707CF |
| MICROMOTOR MC2 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC709CF |
| MICROMOTOR MC3 M40 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC710CF |









MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------------------|---|---|---|---|------------|
| MICROMOTOR MC3 M40 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  <p>Spray regulation Swivel Safety Handpiece</p> <p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p> |  | PC710CFSH |
| MICROMOTOR MC3 M40 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p> |  | PC712CF |
| MICROMOTOR MC3 M40 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  <p>Spray regulation Swivel Safety Handpiece</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p> |  | PC712CFSH |

HIGHSPEED TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------|---|---|--|---|------------|
| HIGHSPEED HSH | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC704CF |
| HIGHSPEED HSH | SKEMA 4 DUO |  |  |  | PC704CFW |
| HIGHSPEED HSH | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC706CF |










HIGHSPEED TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------|--|--|--|--|------------|
| HIGHSPEED HSH | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 |  |   |  | PC704CFSH |
| | STERN S200, 220, 250, 280, 300, 320, 380 | | | | |
| HIGHSPEED HSH | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 |  |   |  | PC706CFSH |
| | STERN S200, 220, 250, 280, 300, 320, 380 | | | | |

SYRINGE TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---|--|-------|--|-------------|------------|
| <p>SYRINGE</p> <p>6F LUZZANI O.F. MINIBRIGHT INOX</p> | <p>ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9</p> <p>STERN S200, 220, 250, 280, 300, 320, 380</p> <p>CASTELLINI SKEMA 5, 6, 8 PUMA ELI</p> | | <p>ITAL GRAY RAL7039</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p> | | PC714EL |
| <p>SYRINGE</p> <p>3F LUZZANI MINILIGHT INOX</p> | <p>ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9</p> <p>STERN S200, 220, 250, 280, 300, 320, 380</p> <p>CASTELLINI SKEMA 5, 6, 8 PUMA ELI</p> | | <p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p> | | PC714C |
| <p>SYRINGE</p> <p>3F LUZZANI MINILIGHT INOX</p> | <p>ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9</p> <p>STERN S200, 220, 250, 280, 300, 320, 380</p> <p>CASTELLINI SKEMA 5, 6, 8 PUMA ELI</p> | | <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p> | | PC714L |










SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|----------------|---|---|--|---|------------|
| SCALER EMS | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC701CF |
| SCALER EMS | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC703CF |
| SCALER PZ5 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC701CP5 |

SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---------------------|---|-------|---------------------------|-------------|------------|
| SCALER PZ5 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI | | | | PC701CP5F |
| SCALER EMS O. F. | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI | | | | PC701CFE |
| SCALER EMS O. F. | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI | | | | PC703CFEL |










SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---|---|---|--|---|------------|
| SCALER PIEZOLIGHT 6 OLD PZ7-SCA3 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC701CF4 |
| SCALER PIEZOLIGHT 6 OLD PZ7-SCA3 | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC703CF4 |
| SCALER SATELEC | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC701CFS |







SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|-------------------------------------|---|-------|--|-------------|------------|
| SCALER SATELEC | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI | | <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 7 mm</p> <p>1440 mm</p> | | PC703CFS |
| SCALER SATELEC NEWTRON LED | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI | | <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 8 mm</p> <p>1250 mm</p> | | PC701CFSL |
| SCALER SATELEC NEWTRON LED | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI | | <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 8 mm</p> <p>1440 mm</p> | | PC703CFSL |

SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--|---|---|--|---|------------|
| SCALER PIEZOSTERIL 5 SURGISON | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC705CPS |
| SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC741C |
| SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  |  | PC741CL |

SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--|---|--|--|--|------------|
| SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 7 mm</p> <p>1440 mm</p> |  | PC742C |
| SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW | ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI |  |  <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 7 mm</p> <p>1440 mm</p> |  | PC742CL |

OPEN END TURBINE DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty

Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N

OPEN END CORDONS TURBINE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--------------------------------|-------------------|-------|---|-------------|------------|
| HIGHSPEED LABORATORY HSH B2 | | | <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>250/1000 mm</p> | | PL042B |
| HIGHSPEED HSH B2/3 | | | <p>Spray regulation Swivel</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p> | | PC334 |
| HIGHSPEED HSH B2/3 | | | <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>200/1000 mm</p> | | PC534 |

OPEN END CORDONS TURBINE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--------------------------|-------------------|-------|---|-------------|------------|
| HIGHSPEED HSH M4 | | | <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p> | | PC444 |
| HIGHSPEED HSH M4 | | | <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>200/1000 mm</p> | | PC544 |
| HIGHSPEED HSH M4 O.F. | | | <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p> | | PC655 |

OPEN END CORDONS TURBINE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---------------------------------|-------------------|-------|--|-------------|------------|
| HIGHSPEED HSH M4 O.F. | | | <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p> | | PC655TK |
| HIGHSPEED HSH CASTELLINI OLD | | | <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>WHITE RAL9010</p> <p>1650 mm - Ø 10,5 mm</p> | | PC306C |
| HIGHSPEED HSH CASTELLINI OLD | | | <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p> | | PC306C/2 |

OPEN END CORDONS TURBINE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---|-------------------|-------|--|-------------|------------|
| HIGHSPEED HSH M4 O.F. | | | <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p> | | PC655TKSH |
| HIGHSPEED HSH M4 O.F. MC3 M40 | | | <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p> | | PC666TM |

OPEN END MICROMOTOR DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty










Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N




BIENAIR MICROMOTORS

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|-------------------------------|-------------------|---|--|---|------------|
| MICROMOTOR MC2 (NO LUX) | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 10 mm</p> |  | PC722 |
| MICROMOTOR MC2 | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 10 mm</p> |  | PC732 |
| MICROMOTOR MC3 | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p> |  | PC666 |










BIENAIR MICROMOTORS

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--------------------------|-------------------|---|--|---|------------|
| MICROMOTOR MC3 HSH | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p> |  | PC666TM |
| MICROMOTOR MC3 | |  |  <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p> |  | PC666SH |
| MICROMOTOR MX2 | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1500 mm - Ø 9,5 mm</p> |  | PC734BXG |

BIENAIR MICROMOTORS

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|-----------------------|----------------------|---|---|---|---------------|
| MICROMOTOR MCX | |  |  <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7038</p> <p>1750 mm - Ø 10,5 mm</p> <p>Spray regulation</p> <p>Swivel</p> |  | PC735MCX |







KAVO MICROMOTORS

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------------------------|-------------------|---|--|---|------------|
| MICROMOTOR K192 | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>2000 mm - Ø 10,5 mm</p> |  | PC124PKN |
| MICROMOTOR K190 K196 K201 | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 10,5 mm</p> |  | PC135PKN |
| MICROMOTOR K700 K701 K703 | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 10,5 mm</p> |  | PC148PKN |




NSK MICROMOTORS

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---|-------------------|-------|--|-------------|------------|
| MICROMOTOR NLX NANO - PLUS ADEC AE53-NLZ | | | <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1800 mm - Ø 8,5 mm</p> | | PC601N |
| MICROMOTOR TIM 40 EL400 | | | <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 9,5 mm</p> | | PC602N |
| MICROMOTOR M40 HSH | | | <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 11 mm</p> | | PC666TM |





NSK MICROMOTORS

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--------------------------|-------------------|---|---|---|------------|
| MICROMOTOR M40 HSH | |  |  <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1700 mm - Ø 11 mm</p> |  | PC901NSH |
| MICROMOTOR NLX-NL400 | |  |  <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 10,5 mm</p> |  | PC603N |







FONA MICROMOTORS

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|----------------------|----------------------|---|--|---|---------------|
| MICROMOTOR BL | |  |  <p>Spray regulation</p> <p>Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 9 mm</p> |  | PC233F |

PLANMECA MICROMOTORS

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---------------------------|----------------------|---|--|---|---------------|
| MICROMOTOR MINETTO | |  | <p>RES. 2,87 KΩ</p> <p>Spray regulation</p> <p>Swivel</p>  <p>https://vimeo.com/480812578</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1800 mm - Ø 10 mm</p>  |  | PC888P |

TEKNEDENTAL MICROMOTORS

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--|-------------------|--|---|---|------------|
| MICROMOTOR TKD DEFINITIVE LED | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 11 mm</p> |  | PC668 |
| MICROMOTOR TKD DEFINITIVE LED | |  |  <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 11 mm</p> |  | PC668SH |

OPEN END SCALER DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

- ⊕ **High quality silicone sheath**
- ⊕ **Latex free**
- ⊕ **12 months warranty**

Technical data:

Maximum pressure supply/return air: 4.5 bar / 65 psi










Maximum pressure air/water spray: 4 bar /58 psi

Tensile strength: 70 N

OPEN END CORDONS SCALER

| HANDPIECE TYPE | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|----------------|-------|---|-------------|------------|
| EMS® | | <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 7 mm</p> | | PC314 |
| EMS® | | <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 7 mm</p> | | PC314/B |
| EMS® O.F. | | <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 8 mm</p> | | PC314L |

OPEN END CORDONS SCALER

| HANDPIECE TYPE | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--|---|--|---|-----------------|
| EMS® I-PIEZON WHITE COLOR |  |  <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 1850 mm - Ø 8 mm</p> |  | PC314IPL |
| SATELEC® NSK VARIOS |  |  <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 7 mm</p> |  | PS115 |
| SATELEC® 3 WIRES |  |  <p>EXTENSION STRING</p> <p>WHITE RAL9010</p> <p>L 2000 mm - Ø 8 mm</p> |  | PS113 |

OPEN END CORDONS SCALER

| HANDPIECE TYPE | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--|-------|---|-----------------------------|------------|
| SATELEC® O.F. ULTRADENT NEWTRON | | <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 8 mm</p> | | PS115L |
| NSK® VARIOS LUX | | <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 1850 mm - Ø 7 mm</p> | <p>with LED 3.5 VDC</p> | PC118N |
| MECTRON® | | <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 1700 mm - Ø 8 mm</p> | | PS117MN |

KAVO DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty

Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N

MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------------------------|--|-------|--|-------------|------------|
| MICROMOTOR K700 K701 K703 | 1057 1060 | | <p>JUMPER 3-4</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> | | PC148KN13 |
| MICROMOTOR K700 K701 K703 | 1058 1062 1063 1065 1066 1078 1080 | | <p>JUMPER 1-3-4</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> | | PC148KN |
| MICROMOTOR K700 K701 K703 | 1058 1062 1063 1065 1066 1078 1080 | | <p>JUMPER 1-3-4</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10,5 mm</p> <p>1450 mm</p> | | PC148KN14 |

MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---|--|-------|---|-------------|------------|
| MICROMOTOR K700 K701 K703 | E30S - E30T E50S - E70S E50TM - E70C E70T - E80C E80T 1058 LIFE | | <p>RES. 4,2 KΩ</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> | | PC148KN5 |
| MICROMOTOR K700 K701 K703 | E30S E30T | | <p>RES. 4,2 KΩ</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> | | PC148KN3 |
| MICROMOTOR K190 K196 K200 TMI MJ660 | 1057 1060 | | <p>JUMPER 3-4</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> | | PC146KN |




MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--|---|-------|---|-------------|-----------------|
| <p>MICROMOTOR</p> <p>K190</p> <p>K196</p> <p>K200</p> <p>TMI MJ660</p> | <p>1058</p> <p>1062</p> <p>1063</p> <p>1065</p> <p>1066</p> <p>1080</p> | | <p>JUMPER 1-3-4</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> | | <p>PC135KN</p> |
| <p>MICROMOTOR</p> <p>K190</p> <p>K196</p> <p>K200</p> <p>TMI MJ660</p> | <p>1066</p> <p>1080</p> | | <p>JUMPER 1-3-4</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1400 mm - Ø 10,5 mm</p> <p>1500 mm</p> | | <p>PC135KNL</p> |

HIGHSPEED TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------|--|-------|--|-------------|------------|
| HIGHSPEED HSH | 1057 1058 1060 1062 1063 1065 1066 1078 1080 | | <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10 mm</p> <p>1350 mm</p> | | PC103KN |
| HIGHSPEED HSH | E50TM E70C E70T E80C E80T 1058 LIFE | | <p>RES. 2,2 KΩ RID SENSOR</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10 mm</p> <p>1350 mm</p> | | PC104KN |
| HIGHSPEED HSH | E50T E70T E80T PENDED | | <p>RES. 2,2 KΩ RID SENSOR</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1550 mm - Ø 10 mm</p> <p>1650 mm</p> | | PC104KN14 |

HIGHSPEED TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|-------------------------------|--|---|---|---|----------------|
| <p>HIGHSPEED HSH</p> | <p>E30S E30T E50S E70S 1058 LIFE S</p> |  | <p>RES. 2,2 KΩ</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10 mm</p> <p>1350 mm</p>  |  | <p>PC105KN</p> |










SYRINGE TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--|---|-------|---|-------------|------------|
| <p>SYRINGE</p> <p>3F KAVO (07730500)</p> | <p>1058</p> <p>1062</p> <p>1063</p> <p>1065</p> <p>1066</p> <p>1080</p> | | <p>LIGHT GRAY RAL7035</p> <p>1150 mm - Ø 10,5 mm</p> <p>1250 mm</p> | | PC168KE |
| <p>SYRINGE</p> <p>4F KAVO (07730000)</p> | <p>1058</p> <p>1062</p> <p>1063</p> <p>1065</p> <p>1066</p> <p>1080</p> | | <p>LIGHT GRAY RAL7035</p> <p>1150 mm - Ø 10,5 mm</p> <p>1250 mm</p> | | PC168KN |
| <p>SYRINGE</p> <p>4F KAVO (07730000)</p> | <p>1066</p> <p>1080</p> | | <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10,5 mm</p> <p>1450 mm</p> | | PC168KNL |

SYRINGE TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---|---|-------|---|-------------|-----------------|
| <p>SYRINGE</p> <p>4F KAVO (07730000)</p> | <p>E50TM E70C E70T E80C E80T 1058 LIFE TM</p> | | <p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> | | <p>PC168KV</p> |
| <p>SYRINGE</p> <p>4F KAVO (07730000)</p> | <p>E50S E70S E70T E80T</p> | | <p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10,5 mm</p> <p>1450 mm</p> | | <p>PC168KVL</p> |
| <p>SYRINGE</p> <p>6F LUZZANI O.F. MINIBRIGHT INOX</p> | <p>E50S E70S E70T E80T</p> | | <p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1495 mm</p> | | <p>PC168KMB</p> |

SYRINGE TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--|--|---|--|---|------------|
| SYRINGE LUZZANI 3F MINILIGHT INOX | 1057 1058 1060 1062 1063 1065 1066 1078 1080 |  |  <p>LIGHT GRAY RAL7035</p> <p>1200 mm - Ø 10,5 mm</p> <p>1445 mm</p> |  | PC179KN |
| SYRINGE LUZZANI 3F MINILIGHT INOX | E50TM E70C E70T E80C E80T 1058 LIFE TM |  |  <p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1495 mm</p> |  | PC179KN2 |
| SYRINGE LUZZANI 3F MINILIGHT INOX | E50S E70S E70T E80T |  |  <p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 10,5 mm</p> <p>1745 mm</p> |  | PC179KN2L |

SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---------------------|---|-------|---|-------------|------------|
| SCALER PIEZOLED | E50S E50TM E70C E70S E70T E80C E80T 1058 LIFE S 1058 LIFE TM 1058LIFE CART | | <p>RES. 1,3 KΩ RID SENSOR</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p> | | PC151KRWG |
| SCALER EMS O. F. | E30S E30T | | <p>RES. 1,3 KΩ</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p> | | PC151KE |
| SCALER EMS O. F. | E30S E30T | | <p>RES. 1,3 KΩ</p> <p>Spray regulation Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1440 mm - Ø 8 mm</p> <p>1520 mm</p> | | PC151KEL |

SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---------------------|---|-------|--|-------------|------------|
| SCALER EMS O. F. | E50S E50TM E70C E70S E70T E80C E80T 1058 LIFE S 1058 LIFE TM 1058LIFE CART | | <p>RES. 1,3 KΩ RID SENSOR</p> <p>Spray regulation </p> <p>Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p> | | PC151KEW |
| SCALER EMS O. F. | 1058 1062 1063 1065 1066 1080 | | <p>JUMPER 1-3</p> <p>Spray regulation </p> <p>Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p> | | PC151KEWG |

PLANMECA DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty

Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N








MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|-------------------|-------------------|-------|--|-------------|------------|
| MICROMOTOR MC3 | | | <p>RES. 83 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>890 mm - Ø 10 mm</p> <p>980 mm</p> | | PC812P |
| MICROMOTOR MC3 | | | <p>RES. 83 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p> | | PC812PL |
| MICROMOTOR MC3 | | | <p>RES. 825 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>890 mm - Ø 10 mm</p> <p>980 mm</p> | | PC812P2 |







MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------------------------|-------------------|-------|---|-------------|------------|
| MICROMOTOR MC3 | | | <p>RES. 825 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p> | | PC812PL2 |
| MICROMOTOR K190 K196 K200 | | | <p>RES. 267 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p> | | PC823PL |
| MICROMOTOR K190 K196 K200 | | | <p>RES. 267 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>900 mm - Ø 10 mm</p> <p>1000 mm</p> | | PC823P |





MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--|-------------------|---|--|---|------------|
| MICROMOTOR K190 K196 K200 | |  |  <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>900 mm - Ø 10 mm</p> <p>1000 mm</p> |  | PC823P/2 |
| MICROMOTOR MINETTO | |  |  <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1800 mm - Ø 10 mm</p> <p></p> <p>https://vimeo.com/480812578</p> |  | PC888P |




HIGHSPEED TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------|-------------------|---|---|---|------------|
| HIGHSPEED HSH | |  | <p>RES. 825 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p> |  | PC801P |
| HIGHSPEED HSH | |  | <p>RES. 430 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>825 mm - Ø 10 mm</p> <p>915 mm</p> |  | PC802P |
| HIGHSPEED HSH | |  | <p>RES. 430 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>840 mm - Ø 10 mm</p> <p>930 mm</p> |  | PC802PL |

HIGHSPEED TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------|-------------------|---|---|---|------------|
| HIGHSPEED HSH | |  | <p>RES. 430 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p> |  | PC803PVL |
| HIGHSPEED HSH | |  | <p>RES. 560 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>915 mm - Ø 10 mm</p> <p>825 mm</p> |  | PC803P |




SYRINGE TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|-------------------------------|----------------------|---|---|---|---------------|
| <p>SYRINGE LUZZANI 3F</p> | |  | <p>RES. 1,78 KΩ</p>  <p>LIGHT GRAY RAL7035</p> <p>910 mm - Ø 10 mm</p> <p>1155 mm</p> |  | <p>PC845P</p> |

SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|----------------------------|-------------------|---|--|---|------------|
| SCALER EMS | |  | <p>RES. 1,33 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>930 mm - Ø 7 mm</p> <p>1020 mm</p> |  | PC834P |
| SCALER EMS O. F. | |  | <p>RES. 1,33 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>930 mm - Ø 8 mm</p> <p>1020 mm</p> |  | PC834PL |
| SCALER SATELEC O. F. | |  | <p>RES. 1,2 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>930 mm - Ø 8 mm</p> <p>1020 mm</p> |  | PC834PSL |

SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|-------------------|----------------------|---|---|---|---------------|
| SCALER SATELEC | |  | <p>RES. 1,2 KΩ</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>930 mm - Ø 7 mm</p> <p>1020 mm</p>  |  | PC834PS |

SIRONA DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty

Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N







MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|--------------------------|--|---|--|---|------------|
| MICROMOTOR MC3 M40 | C6 E3 C3 C3 + |  |  <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1000 mm - Ø 10 mm</p> <p>1100 mm</p> |  | PC233S |
| MICROMOTOR MC3 M40 | C1 C1 + C2 C2 + C4 C4 + C5 + M1 + |  |  <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10 mm</p> <p>1450 mm</p> |  | PC233SL |
| MICROMOTOR EL1 | C1 C1 + C2 C2 + C4 C4 + C5 + M1 + |  |  <p>RES. 4,3 kΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1300 mm - Ø 10 mm</p> <p>1380 mm</p> |  | PC224C4N |










MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|----------------------|------------------------|-------|--|-------------|------------|
| MICROMOTOR EL1 | C3 C3 + C6 | | <p>RES. 4,3 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1000 mm - Ø 10 mm</p> <p>1080 mm</p> | | PC224SN |
| MICROMOTOR EL1 | C8 C8 + PRO FEEL | | <p>RES. 30,4 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10 mm</p> <p>1450 mm</p> | | PC224S/2N |
| MICROMOTOR BL ISO | TENEO | | <p>RES. 2 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 10 mm</p> <p>1370 mm</p> | | PC235SBI |




MICROMOTORS TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------------|---|---|--|---|------------|
| MICROMOTOR BL ISO C | TENEO SINIUS INTEGO INTEGO PRO |  | <p>RES. 3,88 KΩ</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 10 mm</p> <p>1370 mm</p>  |  | PC235SBIC |
| MICROMOTOR BL | TENEO |  | <p>RES. 1,6 KΩ</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 10 mm</p> <p>1370 mm</p>  |  | PC235SBL |







HIGHSPEED TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------|--|---|---|---|------------|
| HIGHSPEED HSH | C6 E3 C3 C3 + |  |  <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1000 mm - Ø 10 mm</p> <p>1080 mm</p> |  | PC204SN |
| HIGHSPEED HSH | A1, A2, A4, A5, AC, C1, C1 +, C2, C2 +, C4, C4 +, C5 +, E1, E2, E4, E5, EC, M1, C8+, PRO FEEL+, TENE0, SINIUS, INTEGO, INTEGO PRO |  |  <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 10 mm</p> <p>1370 mm</p> |  | PC204S/2N |
| HIGHSPEED HSH | C8 PRO FEEL + |  |  <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1300 mm - Ø 10 mm</p> <p>1400 mm</p> |  | PC204S/2E |

HIGHSPEED TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|------------------|--|---|---|---|------------|
| HIGHSPEED HSH | A1, A2, A4, A5, AC, C1, C1 +, C2, C2 +, C4, C4 +, C5 +, E1, E2, E4, E5, EC, M1, M1 96, M1 +, TENEO, SINIUS, INTEGO, INTEGO PRO |  |  <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1560 mm - Ø 10 mm</p> <p>1660 mm</p> |  | PC204S/2L |

SCALER TYPE

| HANDPIECE TYPE | DENTAL UNIT MODEL | IMAGE | LENGHT / DIAMETER / COLOR | CONNECTIONS | ORDER CODE |
|---------------------|---|--|--|--|------------|
| SCALER EMS O. F. | C3 C3+ C6 |  | <p>RES. 47 KΩ</p> <p>Spray regulation Swivel</p>  <p>LIGHT GRAY RAL7035</p> <p>1000 mm - Ø 8 mm</p> <p>1100 mm</p> |  | PC237SEC |
| SCALER EMS O. F. | C8+ PRO FEEL + INTEGO INTEGO PRO |  | <p>Spray regulation Swivel</p>  <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 8 mm</p> <p>1370 mm</p> |  | PC237SEW |