

DENTAL HOSES

 CEFLA	<u>P. 2</u>
 OPEN END TURBINE	<u>P.16</u>
 OPEN END MICROMOTOR	<u>P.21</u>
 OPEN END SCALER	<u>P.31</u>
 KAVO	<u>P.36</u>
 PLANMECA	<u>P.48</u>
 SIRONA	<u>P.57</u>

CEFLA DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty










Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N










MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC2	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC707CF
MICROMOTOR MC2	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC709CF
MICROMOTOR MC3 M40	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC710CF




MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC3 M40	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel Safety Handpiece</p> <p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p>		PC710CFSH
MICROMOTOR MC3 M40	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p>		PC712CF
MICROMOTOR MC3 M40	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel Safety Handpiece</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p>		PC712CFSH










MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR iXR3L iMMrL	ANTHOS A3, 5/+, 6/+, 7+, 9, L6, 9, R7 STERN S200, 250, 280TRc, 380TRc, 300, 220TR, 230TR CASTELLINI SKEMA 5, 6, 8				PC733CX3
MICROMOTOR iXR3L iMMrL	ANTHOS A3, 5/+, 6/+, 7+, 9, L6, 9, R7 STERN S200, 250, 280TRc, 380TRc, 300, 220TR, 230TR CASTELLINI SKEMA 5, 6, 8				PC733CX3L
MICROMOTOR iXS4 iMMs UNTIL 2015	ANTHOS A7+ L9, R7 STERN S320TR, 280TRC, 300, 380TRC CASTELLINI SKEMA 5, 6, 8				PC744CX4







MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>MICROMOTOR</p> <p>iXS4 iMMs</p> <p>UNTIL 2015</p>	<p>ANTHOS A7+ L9, R7</p> <p>STERN S320TR, 280TRC, 300, 380TRC</p> <p>CASTELLINI SKEMA 5, 6, 8</p>				<p>PC744CX4L</p>










HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p>		PC704CF
HIGHSPEED HSH	SKEMA 4 DUO		 <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p>		PC704CFW
HIGHSPEED HSH	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p>		PC706CF










HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9				PC704CFSH
	STERN S200, 220, 250, 280, 300, 320, 380				
HIGHSPEED HSH	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9				PC706CFSH
	STERN S200, 220, 250, 280, 300, 320, 380				










SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>SYRINGE</p> <p>6F LUZZANI O.F. MINIBRIGHT INOX</p>	<p>ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9</p> <p>STERN S200, 220, 250, 280, 300, 320, 380</p> <p>CASTELLINI SKEMA 5, 6, 8 PUMA ELI</p>		<p>ITAL GRAY RAL7039</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p> 		<p>PC714EL</p>
<p>SYRINGE</p> <p>3F LUZZANI MINILIGHT INOX</p>	<p>ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9</p> <p>STERN S200, 220, 250, 280, 300, 320, 380</p> <p>CASTELLINI SKEMA 5, 6, 8 PUMA ELI</p>		<p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 10,5 mm</p> <p>1250 mm</p> 		<p>PC714C</p>
<p>SYRINGE</p> <p>3F LUZZANI MINILIGHT INOX</p>	<p>ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9</p> <p>STERN S200, 220, 250, 280, 300, 320, 380</p> <p>CASTELLINI SKEMA 5, 6, 8 PUMA ELI</p>		<p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 10,5 mm</p> <p>1440 mm</p> 		<p>PC714L</p>










SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER EMS	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CF
SCALER EMS	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC703CF
SCALER PZ5	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CP5

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER PZ5	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CP5F
SCALER EMS O. F.	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CFE
SCALER EMS O. F.	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC703CFEL










SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER PIEZOLIGHT 6 OLD PZ7-SCA3	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CF4
SCALER PIEZOLIGHT 6 OLD PZ7-SCA3	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC703CF4
SCALER SATELEC	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC701CFS







SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER SATELEC	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		<p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 7 mm</p> <p>1440 mm</p>		PC703CFS
SCALER SATELEC NEWTRON LED	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		<p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1160 mm - Ø 8 mm</p> <p>1250 mm</p>		PC701CFSL
SCALER SATELEC NEWTRON LED	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		<p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 8 mm</p> <p>1440 mm</p>		PC703CFSL

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER PIEZOSTERIL 5 SURGISON	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC705CPS
SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC741C
SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI				PC741CL

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 7 mm</p> <p>1440 mm</p>		PC742C
SCALER PIEZOLIGHT 6 PZ7-SCA3 NEW	ANTHOS A3, 5, 6, 7, 9/+ L6, 7, 9 STERN S200, 220, 250, 280, 300, 320, 380 CASTELLINI SKEMA 5, 6, 8 PUMA ELI		 <p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1350 mm - Ø 7 mm</p> <p>1440 mm</p>		PC742CL

OPEN END TURBINE DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty

Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N










OPEN END CORDONS TURBINE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED LABORATORY HSH B2			 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>250/1000 mm</p>		PL042B
HIGHSPEED HSH B2/3			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC334
HIGHSPEED HSH B2/3			 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>200/1000 mm</p>		PC534







OPEN END CORDONS TURBINE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH M4			<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p>		PC444
HIGHSPEED HSH M4			<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>200/1000 mm</p>		PC544
HIGHSPEED HSH M4 O.F.			<p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p>		PC655

OPEN END CORDONS TURBINE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH M4 O.F.			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC655TK
HIGHSPEED HSH CASTELLINI OLD			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>WHITE RAL9010</p> <p>1650 mm - Ø 10,5 mm</p>		PC306C
HIGHSPEED HSH CASTELLINI OLD			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC306C/2

OPEN END CORDONS TURBINE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH M4 O.F.			 <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC655TKSH
HIGHSPEED HSH M4 O.F. MC3 M40			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p>		PC666TM

OPEN END MICROMOTOR DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty










Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N




BIEN AIR MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC2 (NO LUX)			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 10 mm</p>		PC722
MICROMOTOR MC2			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 10 mm</p>		PC732
MICROMOTOR MC3			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC666










BIEN AIR MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC3 HSH			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 11 mm</p>		PC666TM
MICROMOTOR MC3			 <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 10,5 mm</p>		PC666SH
MICROMOTOR MX2			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1500 mm - Ø 9,5 mm</p>		PC734BXG

BIEN AIR MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MCX					PC735MCX










KAVO MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K192			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>2000 mm - Ø 10,5 mm</p>		PC124PKN
MICROMOTOR K190 K196 K201			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 10,5 mm</p>		PC135PKN
MICROMOTOR K700 K701 K703			 <p>Spray regulation Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 10,5 mm</p>		PC148PKN




NSK MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR NLX NANO - PLUS ADEC AE53-NLZ			<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1800 mm - Ø 8,5 mm</p>		PC601N
MICROMOTOR TIM 40 EL400			<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 9,5 mm</p>		PC602N
MICROMOTOR M40 HSH			<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 11 mm</p>		PC666TM




NSK MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR M40 HSH			 <p>Spray regulation Swivel Safety Handpiece</p> <p>EXTENSION STRING</p> <p>ITAL GRAY RAL7038</p> <p>1700 mm - Ø 11 mm</p>		PC901NSH
MICROMOTOR NLX-NL400			 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 10,5 mm</p>		PC603N
MICROMOTOR NBX			 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>2200 mm - Ø 9 mm</p>		PC604N

FONA MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR BL			 <p>Spray regulation</p> <p>Swivel</p> <p>EXTENSION STRING</p> <p>LIGHT GRAY RAL7035</p> <p>1600 mm - Ø 9 mm</p>		PC233F

PLANMECA MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>MICROMOTOR MINETTO</p>					<p>PC888P</p>

TEKNEDENTAL MICROMOTORS

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR TKD DEFINITIVE LED			<p>Spray regulation Swivel</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 11 mm</p>		PC668
MICROMOTOR TKD DEFINITIVE LED			<p>Spray regulation Swivel Safety Handpiece</p> <p>ITAL GRAY RAL7038</p> <p>1650 mm - Ø 11 mm</p>		PC668SH

OPEN END SCALER DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses

- ⊕ **High quality silicone sheath**
- ⊕ **Latex free**
- ⊕ **12 months warranty**

Technical data:

Maximum pressure supply/return air: 4.5 bar / 65 psi

Maximum pressure air/water spray: 4 bar /58 psi

Tensile strength: 70 N










OPEN END CORDONS SCALER

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
EMS®		<p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 7 mm</p>		PC314
EMS®		<p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 7 mm</p>		PC314/B
EMS® O.F.		<p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 8 mm</p>		PC314L




OPEN END CORDONS SCALER

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
EMS® O.F.		<p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 1980 mm - Ø 8 mm</p>		PC314LW
EMS® I-PIEZON WHITE COLOR		<p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 1850 mm - Ø 8 mm</p>		PC314IPL
SATELEC® NSK VARIOS		<p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 7 mm</p>		PS115

OPEN END CORDONS SCALER

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SATELEC® 3 WIRES		 <p>EXTENSION STRING</p> <p>WHITE RAL9010</p> <p>L 2000 mm - Ø 8 mm</p>		PS113
SATELEC® O.F. ULTRADENT NEWTRON		 <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 2000 mm - Ø 8 mm</p>		PS115L
NSK® VARIOS LUX		 <p>EXTENSION STRING</p> <p>LIGHT GREY RAL7035</p> <p>L 1850 mm - Ø 7 mm</p>	 <p>with LED 3.5 VDC</p>	PC118N

OPEN END CORDONS SCALER

HANDPIECE TYPE	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MECTRON®				PS117MN

KAVO DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty



















Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N

MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR K700 K701 K703	1057 1060		<p>JUMPER 3-4</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> 		PC148KN13
MICROMOTOR K700 K701 K703	1058 1062 1063 1065 1066 1078 1080		<p>JUMPER 1-3-4</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> 		PC148KN
MICROMOTOR K700 K701 K703	1058 1062 1063 1065 1066 1078 1080		<p>JUMPER 1-3-4</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10,5 mm</p> <p>1450 mm</p> 		PC148KN14

MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>MICROMOTOR</p> <p>K700 K701 K703</p>	<p>E30S - E30T E50S - E70S E50TM - E70C E70T - E80C E80T 1058 LIFE</p>		<p>RES. 4,2 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> 	 	PC148KN5
<p>MICROMOTOR</p> <p>K700 K701 K703</p>	<p>E30S E30T</p>		<p>RES. 4,2 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> 	 	PC148KN3
<p>MICROMOTOR</p> <p>K190 K196 K200 TMI MJ660</p>	<p>1057 1060</p>		<p>JUMPER 3-4</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p> 	 	PC146KN




MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>MICROMOTOR</p> <p>K190 K196 K200 TMI MJ660</p>	<p>1058 1062 1063 1065 1066 1080</p>		<p>JUMPER 1-3-4</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p>		<p>PC135KN</p>
<p>MICROMOTOR</p> <p>K190 K196 K200 TMI MJ660</p>	<p>1066 1080</p>		<p>JUMPER 1-3-4</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1400 mm - Ø 10,5 mm</p> <p>1500 mm</p>		<p>PC135KNL</p>










HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH	1057 1058 1060 1062 1063 1065 1066 1078 1080		<p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10 mm</p> <p>1350 mm</p>		PC103KN
HIGHSPEED HSH	E50TM E70C E70T E80C E80T 1058 LIFE		<p>RES. 2,2 KΩ RID SENSOR</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10 mm</p> <p>1350 mm</p>		PC104KN
HIGHSPEED HSH	E50T E70T E80T PENDED		<p>RES. 2,2 KΩ RID SENSOR</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1550 mm - Ø 10 mm</p> <p>1650 mm</p>		PC104KN14










HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>HIGHSPEED HSH</p>	<p>E30S E30T E50S E70S 1058 LIFE S</p>		 <p>RES. 2,2 KΩ</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10 mm</p> <p>1350 mm</p>		<p>PC105KN</p>










SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>SYRINGE</p> <p>3F KAVO (07730500)</p>	<p>1058</p> <p>1062</p> <p>1063</p> <p>1065</p> <p>1066</p> <p>1080</p>				PC168KE
<p>SYRINGE</p> <p>4F KAVO (07730000)</p>	<p>1058</p> <p>1062</p> <p>1063</p> <p>1065</p> <p>1066</p> <p>1080</p>				PC168KN
<p>SYRINGE</p> <p>4F KAVO (07730000)</p>	<p>1066</p> <p>1080</p>				PC168KNL
















SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>SYRINGE</p> <p>4F KAVO (07730000)</p>	<p>E50TM E70C E70T E80C E80T 1058 LIFE TM</p>		<p>RES. 750 Ω</p>  <p>1250 mm - Ø 10,5 mm</p> <p>1350 mm</p>		<p>PC168KV</p>
<p>SYRINGE</p> <p>4F KAVO (07730000)</p>	<p>E50S E70S E70T E80T</p>		<p>RES. 750 Ω</p>  <p>1350 mm - Ø 10,5 mm</p> <p>1450 mm</p>		<p>PC168KVL</p>
<p>SYRINGE</p> <p>6F LUZZANI O.F. MINIBRIGHT INOX</p>	<p>E50S E70S E70T E80T</p>		<p>RES. 750 Ω</p>  <p>1250 mm - Ø 10,5 mm</p> <p>1495 mm</p>		<p>PC168KMB</p>

SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>SYRINGE</p> <p>LUZZANI 3F MINILIGHT INOX</p>	<p>1057 1058 1060 1062 1063 1065 1066 1078 1080</p>		 <p>LIGHT GRAY RAL7035</p> <p>1200 mm - Ø 10,5 mm</p> <p>1445 mm</p>		<p>PC179KN</p>
<p>SYRINGE</p> <p>LUZZANI 3F MINILIGHT INOX</p>	<p>E50TM E70C E70T E80C E80T 1058 LIFE TM</p>		 <p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 10,5 mm</p> <p>1495 mm</p>		<p>PC179KN2</p>
<p>SYRINGE</p> <p>LUZZANI 3F MINILIGHT INOX</p>	<p>E50S E70S E70T E80T</p>		 <p>RES. 750 Ω</p> <p>LIGHT GRAY RAL7035</p> <p>1500 mm - Ø 10,5 mm</p> <p>1745 mm</p>		<p>PC179KN2L</p>

SCALER TYPE


HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER* PIEZOLED EMS O. F.	E50S E50TM E70C E70S E70T E80C E80T 1058 LIFE S 1058 LIFE TM 1058LIFE CART		<p>RES. 1,3 KΩ RID SENSOR</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p> 		PC151KRWG
SCALER EMS O. F.	E30S E30T		<p>RES. 1,3 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p> 		PC151KE
SCALER EMS O. F.	E30S E30T		<p>RES. 1,3 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1440 mm - Ø 8 mm</p> <p>1520 mm</p> 		PC151KEL

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER EMS O. F.	E50S E50TM E70C E70S E70T E80C E80T 1058 LIFE S 1058 LIFE TM 1058LIFE CART		<p>RES. 1,3 KΩ RID SENSOR</p> <p>Spray regulation </p> <p>Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p>		PC151KEW
SCALER EMS O. F.	1058 1062 1063 1065 1066 1080		<p>JUMPER 1-3</p> <p>Spray regulation </p> <p>Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p>		PC151KEWG

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER * EMS O. F. KAVO PIEZOLED	E50S E50TM E70C E70S E70T E80C E80T 1058		<p>RES. 1,3 KΩ</p> <p>JUMPER 1-3</p> <p>Spray regulation </p> <p>Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1255 mm - Ø 8 mm</p> <p>1350 mm</p>	<p>REVERSE POLARITY</p>	PC152KEWG
				PC152KEWGI	
SCALER * EMS O. F. KAVO PIEZOLED	E30S E30T		<p>RES. 1,3 KΩ</p> <p>Spray regulation </p> <p>Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>1250 mm - Ø 8 mm</p> <p>1350 mm</p>		PC153KRWG



PLANMECA DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses







Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty










Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N







MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC3			<p>RES. 83 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>890 mm - Ø 10 mm</p> <p>980 mm</p>		PC812P
MICROMOTOR MC3			<p>RES. 83 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p>		PC812PL
MICROMOTOR MC3			<p>RES. 825 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>890 mm - Ø 10 mm</p> <p>980 mm</p>		PC812P2







MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC3			<p>RES. 825 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p> 		PC812PL2
MICROMOTOR K190 K196 K200			<p>RES. 267 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p> 		PC823PL
MICROMOTOR K190 K196 K200			<p>RES. 267 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>900 mm - Ø 10 mm</p> <p>1000 mm</p> 		PC823P





MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>MICROMOTOR</p> <p>K190 K196 K200</p>					<p>PC823P/2</p>
<p>MICROMOTOR</p> <p>MINETTO</p>					<p>PC888P</p>




HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH			<p>RES. 825 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p>		PC801P
HIGHSPEED HSH			<p>RES. 430 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>825 mm - Ø 10 mm</p> <p>915 mm</p>		PC802P
HIGHSPEED HSH			<p>RES. 430 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>840 mm - Ø 10 mm</p> <p>930 mm</p>		PC802PL













HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH			<p>RES. 430 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>925 mm - Ø 10 mm</p> <p>1015 mm</p>		PC803PVL
HIGHSPEED HSH			<p>RES. 560 Ω</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>915 mm - Ø 10 mm</p> <p>825 mm</p>		PC803P




SYRINGE TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>SYRINGE LUZZANI 3F</p>			<p>RES. 1,78 KΩ</p>  <p>LIGHT GRAY RAL7035</p> <p>910 mm - Ø 10 mm</p> <p>1155 mm</p>		<p>PC845P</p>

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER EMS			<p>RES. 1,33 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>930 mm - Ø 7 mm</p> <p>1020 mm</p>		PC834P
SCALER EMS O. F.			<p>RES. 1,33 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>930 mm - Ø 8 mm</p> <p>1020 mm</p>		PC834PL
SCALER SATELEC O. F.			<p>RES. 1,2 KΩ</p> <p>Spray regulation  Swivel </p> <p>LIGHT GRAY RAL7035</p> <p>930 mm - Ø 8 mm</p> <p>1020 mm</p>		PC834PSL

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER SATELEC			<p>RES. 1,2 KΩ</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>930 mm - Ø 7 mm</p> <p>1020 mm</p> 		PC834PS

SIRONA DENTAL HOSES

The dental hoses for Tecnomed Italia dental units are made of silicone; they are Latex Free and comply with the UNI CEI EN ISO 13485 quality system.

The chemically stable gold-plated contacts ensure high electrical conductivity and maximum handpiece performance.

The supply tubes comply with REGULATION (EU) No. 10/2011 concerning plastic materials objects intended to meet food products and the exhaust tubes are resistant to the lubrication oils.

Tecnomed Italia dental hoses meet the strictest medical standards, complying with the essential requirements of the European Directives (93/42 / EEC) 2007/47 / EC on medical devices and 2000/95 / EC RoHs on the restriction of the use of certain specific substances in electrical and electronic equipment.

A complete range of dental hoses:

- Turbines
- Micromotors
- Scalers
- Syringes
- Freely Adaptable Dental Hoses










Features

- ⊕ High quality silicone sheath.
- ⊕ Latex free.
- ⊕ 12 months warranty







Technical data

- ⊕ Maximum pressure supply/return air 4,5 bar / 65 psi
- ⊕ Maximum pressure air/water spray 4 bar / 58 psi
- ⊕ Tensile strenght 70N









MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR MC3 M40	C6 E3 C3 C3 +		 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1000 mm - Ø 10 mm</p> <p>1100 mm</p>		PC233S
MICROMOTOR MC3 M40	C1 C1 + C2 C2 + C4 C4 + C5 + M1 +		 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10 mm</p> <p>1450 mm</p>		PC233SL
MICROMOTOR EL1	C1 C1 + C2 C2 + C4 C4 + C5 + M1 +		 <p>RES. 4,3 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1300 mm - Ø 10 mm</p> <p>1380 mm</p>		PC224C4N










MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR EL1	C3 C3 + C6		<p>RES. 4,3 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1000 mm - Ø 10 mm</p> <p>1080 mm</p>		PC224SN
MICROMOTOR EL1	C8 C8 + PRO FEEL		<p>RES. 30,4 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1350 mm - Ø 10 mm</p> <p>1450 mm</p>		PC224S/2N
MICROMOTOR BL ISO	TENEO		<p>RES. 2 KΩ</p> <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 10 mm</p> <p>1370 mm</p>		PC235SBI




MICROMOTORS TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
MICROMOTOR BL ISO C	TENEO SINIUS INTEGO INTEGO PRO		<p>RES. 3,88 KΩ</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 10 mm</p> <p>1370 mm</p> 	 	PC235SBIC
MICROMOTOR BL	TENEO		<p>RES. 1,6 KΩ</p> <p>Spray regulation</p> <p>Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 10 mm</p> <p>1370 mm</p> 	 	PC235SBL







HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
HIGHSPEED HSH	C6 E3 C3 C3 +		 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1000 mm - Ø 10 mm</p> <p>1080 mm</p>		PC204SN
HIGHSPEED HSH	A1, A2, A4, A5, AC, C1, C1 +, C2, C2 +, C4, C4 +, C5 +, E1, E2, E4, E5, EC, M1, C8 +, PRO FEEL +, TENEQ, SINIUS, INTEGO, INTEGO PRO		 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 10 mm</p> <p>1370 mm</p>		PC204S/2N
HIGHSPEED HSH	C8 PRO FEEL +		 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1300 mm - Ø 10 mm</p> <p>1400 mm</p>		PC204S/2E

HIGHSPEED TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGTH / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
<p>HIGHSPEED HSH</p>	<p>A1, A2, A4, A5, AC, C1, C1 +, C2, C2 +, C4, C4 +, C5 +, E1, E2, E4, E5, EC, M1, M1 96, M1 +, TENEO, SINIUS, INTEGO, INTEGO PRO</p>		 <p>Spray regulation Swivel</p> <p>LIGHT GRAY RAL7035</p> <p>1560 mm - Ø 10 mm</p> <p>1660 mm</p>		<p>PC204S/2L</p>

SCALER TYPE

HANDPIECE TYPE	DENTAL UNIT MODEL	IMAGE	LENGHT / DIAMETER / COLOR	CONNECTIONS	ORDER CODE
SCALER EMS O. F.	C3 C3+ C6		<p>RES. 47 KΩ</p> <p>Spray regulation Swivel</p>  <p>LIGHT GRAY RAL7035</p> <p>1000 mm - Ø 8 mm</p> <p>1100 mm</p>		PC237SEC
SCALER EMS O. F.	C8+ PRO FEEL + INTEGO INTEGO PRO		<p>Spray regulation Swivel</p>  <p>LIGHT GRAY RAL7035</p> <p>1270 mm - Ø 8 mm</p> <p>1370 mm</p>		PC237SENW